

Protein Characterization of Dehydrated Placenta Membrane Products: Artacent Wound and Artacent AC



Babak Safavieh, PhD, Olivia Logan, BS, Denae Landry, BS, Isabella Sledge, MD, Mora Melican, PhD

This poster was first presented at the Symposium on Advanced Wound Care (SAWC Spring 2023) in National Harbor, MD.

Introduction

Human Amniotic Membrane (hAM) is a readily available tissue that is widely used as a wound covering in both chronic and acute wounds. Dehydration is commonly used as a processing step to create a product that is stable at room temperature and has a long shelf life. Two dehydrated membrane products that are widely being used in the clinic are: (i) Artacent Wound®, that consists of bi-layer of amnion dried membrane, and (ii) Artacent AC®, which is tri-layer composed of chorion held between two layers of amnion. In general, dehydrated hAMs have been reported [1] to contain a variety of proteins, growth factors and cytokines which help to not only promote and stimulate cell proliferation and migration, but also alter gene and protein expression profiles of cells. Here, we characterize the protein content in Artacent Wound® and Artacent AC® using Microarray. The statistical analysis will be conducted to determine the combination of biomarkers that are present in various products, and to compare the biomarkers in Artacent products to that of clinical controls. In addition, immunohistochemistry (IHC) was used to demonstrate various protein composition of extracellular matrix (ECM), growth factors, and immune modulating cytokines. Finally, cell attachment and proliferation of human fibroblasts and keratinocytes was assessed in vitro on Artacent AC® and Artacent Wound® products.

Materials and Methods

- 1


HISTOLOGY
A 20 x 20 mm section of each tissue was prepared and sent to Histowiz Inc. for histological analysis. Samples were paraffin embedded and stained with Hematoxylin and Eosin (H&E) and immunohistochemistry (IHC) labels of Collagen I, Collagen III, Collagen V, Collagen VI, Fibronectin, and VEGF for different extra cellular matrix proteins and chemokine cues.
- 2

MICROARRAY
Quantibody® Human Cytokine Antibody Array 1000 assay from RayBiotech Inc. The assay is the multiplexed ELISA quantitative platform enabling to detect and quantify up to 1000 human cytokines in a single experiment. Briefly, samples were incubated on captured antibody for 2 hours. The biotinylated antibody cocktail were added and then incubated again for another 2 hours following by streptavidin labeled dye. The result were then analyzed using fluorescent microarray reader.
- 3

CELL CULTURE
Fibroblast L929 cells were cultured and 10⁴5 cells/ml were added to each sample of each well with 3 mm diameter. Acti-stain™ 488 phalloidin were added and fixed before microscopy imaging as described previously [2].
- 4

STATISTICAL ANALYSIS
Statistical analysis on all results was performed using Two-way ANOVA with the test for multiple comparison using Tukey statistical hypothesis testing group of products of all biomarkers. Statistical significance threshold was set at 0.05 (p < 0.05). Error bars represent a standard deviation of the mean (n = 5). An asterisk shows statistically significant magnitude change compared to all other datasets.

Results



Discussion

Protein microarray analysis and quantification shows Artacent AC and Artacent Wound have various protein, biomarker, and chemical cues for angiogenic and immune modulating cytokines and regenerative biomarkers (Figure 1A). There are no significant differences in angiogenic and immunomodulating biomarkers in Artacent AC, Artacent Wound, and clinical control (Figure 1B&C). These products posed statistical differences in BMP-5, IGFBP-3 and IGFBP-4. Artacent AC has more BMP-5 compared to Artacent Wound and clinical control (P<0.01). BMP-5 is beneficial to promote wound healing and tissue regeneration [3].

Artacent Wound protein composition posed uniform distribution of immune, angiogenic, and regenerative molecules. However, the protein distribution in Artacent AC showed more regenerative growth factors compared to angiogenic and immune modulating cytokines.

Histology H&E and IHC of various ECM biomarkers showed that Artacent Wound consists of various collagen types (I, III, V and VI), in addition to fibronectin and VEGF growth factor. This ECM composition makes the product unique to promote re-epithelization in the wound healing process. Biocompatibility assay of Artacent products improves fibroblast proliferation in comparison with clinical control samples. This may attribute to the fact that Artacent products consist of the enriched biomarker and cytokine molecules.

Conclusions

- Artacent AC shows higher level of BMP-5 biomarker compared to control sample.
- Histological analysis showed Artacent Wound has various ECM proteins and growth factor:
 - Collagen I, III, V, VI
 - Fibronectin
 - VEGF
- Both Artacent AC and Artacent Wound showed better cell proliferation within 24 hours when compared to a similar product as well as collagen sheet.

REFERENCES

- McQuilling JP et al. Int Wound J. 2019;16:827-840.
- Acti-Stain™ Fluorescent Phalloidins. <https://www.cytoskeleton.com/actin/acti-stain>
- Carmagnola D, Pellegrini G, Dellavia C, Rimondini L, Varoni E. Tissue engineering in periodontology: Biological mediators for periodontal regeneration. The International Journal of Artificial Organs. 2019;42(5):241-257. doi:10.1177/0391398119828558

ACKNOWLEDGEMENT

This poster work was supported financially by Tides Medical Research and Development sources. The authors also extend their gratitude to Andr e Gonsoulin and Sarah Marr from Tides Medical and Mitchell Sanders, Mia Hanna, Nikole Siegmund, Lindsay Poland and Da Wei from Alira Health. Additionally, the authors thank prof. Harini Sundararaghavan from Wayne State University.