


# Characterization of Human iPSC-derived Induced Excitatory Neurons for Disease Modeling and Drug Discovery



Katherine Czysz, Robert Bradley, Junyi Ma, Ali Fathi, Rebecca Fiene, Clifford Hogan, Coby Carlson, Simon Hilcove, Scott Schachtele, and Jing Liu


FUJIFILM Cellular Dynamics, Inc., Madison, WI USA




## Abstract

Human induced pluripotent stem cell (iPSC)-derived neurons are integral for elucidating mechanisms and therapeutic targets underlying neurodevelopmental and neurodegenerative disease. Directed differentiation protocols have largely been employed for generating iPSC-derived neural models, because of the ability to generate specific cortical glutamatergic neurons for use in human drug discovery and screening. While these cells are highly predictive for drug screening, directed differentiation methods are challenged to meet the economies of scale required for high-throughput compound testing within preclinical phases of the drug development pipeline. Induction of neural differentiation using NGN2 forward reprogramming of iPSCs offers a robust method for generating scalable quantities of excitatory neurons with low lot-to-lot variability. In this study we utilize NGN2 overexpression, under a doxycycline (DOX) promoter, to forward program highly pure excitatory glutamatergic neurons from iPSC lines with an apparently healthy normal (AHN) background or heterozygous pathogenic R493X nonsense mutation in the granulin gene (GRN) to model frontotemporal dementia (FTD). We first demonstrate that induced excitatory AHN and GRN (R493X) neurons display gene expression profiles similar to published NGN2 protocols. We further characterized these cells, showing the presence of highly pure neuronal populations (> 90% TUJ1-positive) that can be recovered from cryopreservation to form consistent neural cultures without continued use of DOX. When co-cultured with isogenic iPSC-derived astrocytes, both AHN and GRN (493X) induced excitatory neurons produced robust neural networks on multielectrode array platforms. Last, induced excitatory neurons were co-cultured with iPSC-derived astrocytes on ultra-low adhesion 96-well round-bottom plates to form neurospheres, which developed dynamic calcium oscillations as early as 14 days in culture. These data demonstrate a robust and scalable method using NGN2 forward programming to consistently produce fully differentiated induced excitatory neurons. These cells function across a variety of applications that support high-throughput preclinical drug discovery and compound screening.


## Development and Characterization of iPSC-derived Induced Neurons




## Functional Characterization of Induced Excitatory Neurons on MEA




## Analysis of Induced Excitatory Neurons on High-density MEA



## Calcium Waveform Oscillations from Induced Neuron 3D Neurospheres



## Disease Modeling with Induced Neurons: GRN R493X HZ KO



## Summary

These data demonstrate a highly consistent and easy-to-use induced (“forward programmed”) neuron that allows for high-throughput analysis of toxicity and neuronal function. iCell Induced Excitatory Neurons will function without additional differentiation by the user. We also show these cells can be combined with CRISPR engineering to create robust disease models within glutamatergic neurons.

