

Automated In Vitro MDCKII-BCRP Monolayer Assay with Ultra-Fast LC-MS/MS

E. Tímár¹, V. Juhász¹, Z. Urban¹, Zs. Nerada¹, O. Csíkvári¹, Z. Timár¹, A. Németh¹, P. Szoleczky¹

¹SOLVO Biotechnology – a Charles River Company

1 INTRODUCTION

Bidirectional permeability assay:

- In vitro modeling of pharmacokinetic barriers (e.g., Intestinal epithelium, blood-brain-barrier, etc.).
- Informs us about the compounds of interest active and passive transport characteristics.
- Tight cell monolayer grown on a porous support
- MDCKII (Madin-Darby canine kidney strain II) cells overexpressing the BCRP (Breast Cancer Resistant Protein, ABCG2) efflux transporter in their apical membrane (MDCKII-BCRP) and acts as a B (basolateral) to A (apical) pump for its substrates.




Fig 1. In vitro bidirectional permeability assay design for assessing apparent permeability (P_{app}) of a compound.

Challenges of a high throughput system:

- Changing from the typically used 24-well to 96 well plate format.
- High requirements for analytical instrumentation. As sample complexity increases, sensitivity and reliability of MS-based plate reading approaches decrease.

2 METHODS AND MATERIALS

Assay Setup

- Cell seeding was performed 4 days before the experiment onto Millicell-96 culture plates.
- The reagent mix containing the test articles and controls were prepared on a deep well plate (helper plate).
- Proper cell function was validated by analyzing A-B and B-A movement of radiolabeled prazosin, a known BCRP substrate (accepted if Efflux Ratio (ER) >8).
- Cell monolayer integrity was confirmed via measurements with Lucifer yellow (LY) as low permeability control (accepted if $P_{app} < 2 \times 10^{-6}$ cm/s) in each well and metoprolol as high permeability control (n=3).




Fig 2. Bidirectional permeability assay plate layout. The apical and basolateral wells of the assay plates were filled with the reaction mixture prepared on the helper plate. After the proper incubation time the samples were collected on separate microplates.

3 WORKFLOW AND SYSTEM DESIGN




Fig 3. Steps of the experimental workflow and the related instruments.

The assay was performed in a 96-well format plate. The cell seeding and media exchanges were performed manually.

The helper plate was prepared using a Tecan D300 digital dispenser and a Hamilton MicroLab STAR automated liquid handler. The liquid handler, equipped with a 96 CO-RE Head was used for the compound dosing and sample collection. The instrument is also integrated with a BioTek 405 TS microplate washer for the buffer exchange, and a Thermo Cytomat 2 C-LIN incubator.




Fig 4. Schematic representation of the Hamilton MicroLab Star deck layout

LC-MS/MS analysis was performed on an ADDA (Apricot Designs Dual Arm) autosampler coupled to a Sciex 6500+ triple quadrupole mass spectrometer. Method development and sample analysis were performed using LeadScape software, designed for high throughput analysis.

4 RESULTS

	Manual workflow			Automated workflow		
	P_{app} , A-B ($\times 10^6$ cm/s)	P_{app} , B-A ($\times 10^6$ cm/s)	ER	P_{app} , A-B ($\times 10^6$ cm/s)	P_{app} , B-A ($\times 10^6$ cm/s)	ER
	Average \pm SEM	Average \pm SEM	Average \pm EP	Average \pm SEM	Average \pm SEM	Average \pm EP
Prazosin	2.07 \pm 0.29	67.78 \pm 12.62	32.77 \pm 7.59	1.77 \pm 0.17	114.84 \pm 5.52	64.77 \pm 5.52
Lucifer Yellow	1.28 \pm 1.97	-	-	0.237 \pm 0.018	0.370 \pm 0	-
Metoprolol	39.66 \pm 24.50	-	-	41.99 \pm 2.58	59.28 \pm 0.47	1.41 \pm 2.62

Table 1. Bidirectional permeability of control compounds using the manual or automated workflow, respectively.




Fig 5. Efflux ratios (ER) calculated for all compounds tested using the automated assay workflow.

5 CONCLUSIONS

Our results demonstrate that the presented automated monolayer assay system can yield bidirectional permeability data that reliably reproduces values obtained in the manually conducted assay workflow. Permeability of both the high- and low-permeability reference compounds as well as the selected BCRP probe substrate was similar using the two workflows and can serve as appropriate controls when testing unknown substances. The application of an automated assay system where possible can reduce human error, increase throughput and save time, whilst producing reliable and reproducible data.

[1] Agnès Poirier, Renée Portmann, Anne-Christine Cascais, Urs Bader, Isabelle Walter, Mohammed Ullah, and Christoph Funk. The Need for Human Breast Cancer Resistance Protein Substrate and Inhibition Evaluation in Drug Discovery and Development: Why, When, and How? Drug Metab Dispos 42:1466–1477, September 2014
 [2] Houfu Liu, Liang Huang, Yi Li, Tingting Fu, Xueying Sun, Yanyan Zhang, Ruina Gao, Qingfang Chen, Wandong Zhang, Jasinder Sahi, Scott Summerfield, and Kelly Dong. Correlation Between Membrane Protein Expression Levels and Transcellular Transport Activity for Breast Cancer Resistance Protein. DMD Fast Forward. Published on February 16, 2017
 [3] Lukas Cerveny, Petr Pavek, Jana Malakova, Frantisek Staud, and Zdenek Fendrich. Lack of Interactions between Breast Cancer Resistance Protein (BCRP/ABCG2) and Selected Antiepileptic Agents Epilepsia. 47(3):461–468, 2006

