

# High content imaging is now high throughput: End-to-end high content imaging and analysis of autophagy on a subcellular level in minutes



Matthew Boisvert, Will Marshall, Josh Kieler, Amanda Clement

## Measuring autophagic flux at high throughput speeds

- Cells eliminate intracellular waste and defective components through autophagy ("self-eating")
- A key therapeutic target for aging-related dysfunction, autophagy is implicated in neurodegeneration and cancer
- LC3b aggregation measures autophagic flux, resulting from increases in autophagy (stress) or lysosomal blockage
- High content imaging (HCI) allows for cell-level resolution, visualizing then counting aggregates within the cell
- Resolving aggregated LC3b is essential, cell lines can express endogenous LC3b in non-aggregated form
- Screening for autophagic flux requires speed for a many compound library, low variability and accurate detection
- This autophagy assay with the Araceli Platform shows:
  - Submicron resolution to accurately visualize aggregates
  - Elimination of variability with whole well imaging
  - <10-minute scan times for 96, 384 and 1536 plates
  - Comparable analysis to CellProfiler in <1/10 the time

**ARACELI Voyager**  
cell image analysis software

Grow cells: 48 hrs  
Treat cells: 12 hrs  
Autophagy staining  
Fix and stain: 4 hrs  
Set up: <5 minutes  
Imaging: 7-10 min/plate  
Compare analysis pipelines: 30min-3 days

## Methods


### Cell culture and autophagic flux assay

Human bone carcinoma (U2OS) cells were plated at 12500 cells/well in 96 well Grenier CellView plates and cultured for 48 hours in standard conditions (complete media, 5% CO<sub>2</sub>, 37°C). Cells were treated with compound or DMSO in a 1:2 dilution series, 3 replicates/plate, for 12-16 hours. Cells were then fixed in 10% formalin for 15 minutes, permeabilized with 0.1% Triton x-100, then stained and counterstained for 60min/each in 1% BSA, with 1:1000 Rb anti-LC3b (L10382) and 1:1000 Gt anti-Rb Alexa 488 with 1:750 Hoechst and 1:1000 iFluor Phalloidin-647. Cells were washed 4x between steps and all solutions are in PBS pH 7.4.

Treatments: **DMSO**: vehicle control; **[hydroxy]chloroquine**: autophagosome degradation inhibitor; **verapamil**: Ca<sup>2+</sup> inhibitor; **MG-132**: inhibit proteasome

### Comparative analyses


Cytoplasmic spot detection and counting done in parallel between CellProfiler and Araceli Voyager™. For both, nuclei and cell segmented (see figure for Voyager details and example), with nuclei subtracted from cell to define cytoplasm, and spots counted within cytoplasm. For CellProfiler, images were preprocessed to eliminate illumination profile and background intensity offset and spots enhances using difference of gaussians; Voyager automatically flatfield/vignette corrects image from Endeavor. CellProfiler optimization of parameters done on test set of 3 images with full pipeline generating new image files, which were then visually assessed and rerun with new parameters. Voyager's assay development mode used to visualize and tweak each parameter stepwise with near instant analysis feedback. Z' calculated from effect size peak, compared to DMSO control.




### Full well, full plate high resolution scans in <10 minutes

96 well plates imaged at 0.27µm/pixel, 3 channels, with (A) full wells (~85% of well area) collected. Exposures set up to maximize dynamic range at drug concentration peak for each plate, 10µM chloroquine shown with Hoechst (nuclei, B), LC3B (autophagosomes, C), and phalloidin (actin, D); merged image (E).

## Autophagic flux consistently detected with comparable levels in Voyager and CellProfiler



Z' Score	Chloroquine	Hydrochloroquine	MG-132	Verapamil
Voyager	0.8314	0.8503	0.8342	0.6325
CellProfiler	0.9037	0.8355	0.8799	0.9489




## Voyager analysis workflow saves days of assay set up, hours of compute time per plate

Analysis	Set up time	Analysis time	CPU/RAM (gb)
CellProfiler	2 days	4-5 hrs / plate	128 / 256
Voyager	< 15 min	10-20 min / plate	32 / 64

### Voyager workflow

Crop of 5µM chloroquine treatment in Voyager assay development mode:

- Step 1:** Plates loaded into Voyager with automatic flat field/vignette correction;
- Step 2:** Identify and segment nuclei, using machine vision on Hoechst stain;
- Step 3:** Segment cells using actin stain, subtracting nuclei to create cytoplasmic mask;
- Step 4:** Identify spots with template match, optionally segment spots if area and intensity desired.



## HCS and analysis at a submicron subcellular resolution at high throughput speeds

- Autophagic flux reliably detected with several known modulators in dilution
- Good Z' indicates robust effect size, minimal variability
- Plate to results delivered in under 30 minutes
- Araceli Voyager yield equivalent analysis to CellProfiler with:
  - 10% compute time at 25% compute power
  - Easy, straightforward and fast set up
  - Near-instant feedback to refine each step, and assay overall
- 0.27µm pixel size allows for spot detection at submicron scale
- Assay broadly generalizable to other quantitative spot-based HCS assays such as: subcellular protein localization, plaque formation, aggregation, FISH, phagocytosis

## REFERENCES

Djadjadikerta A, Keshri S, Pavel M, Prestil R, Ryan L, Rubinsztein DC. Autophagy Induction as a Therapeutic Strategy for Neurodegenerative Diseases. J Mol Biol. 2020 Apr 3;432(8):2799-2821. doi: 10.1016/j.jmb.2019.12.035.  
 Joachim J, Jiang M, McKnight NC, Howell M, Toose SA. High-throughput screening approaches to identify regulators of mammalian autophagy. Methods. 2015 Mar;75:96-104. doi: 10.1016/j.jmeth.2015.02.002  
 Sharif MN, Mowers EE, Drake LE, Macleod KF. Measuring autophagy in stressed cells. Methods Mol Biol. 2015;1292:129-50. doi: 10.1007/978-1-4939-2522-3\_10.  
 Marinković M, Sprung M, Bujubašić M, Novak I. Autophagy Modulation in Cancer: Current Knowledge on Action and Therapy. Oxid Med Cell Longev. 2018 Jan 31;2018:8023821. doi: 10.1155/2018/8023821