

## Introduction


Biofilms are the leading cause in the stalling of wound healing in clinical settings as they can penetrate deep into tissue and evade topical antimicrobial treatment. *Ex vivo* human skin models allow for moderately high throughput screening while remaining clinically relevant by using a host-derived substrate, allowing for multiple experimental outputs and data types. **This study aimed to determine the biofilm disruption and antibiofilm depth penetration properties of two new prototype wound dressings, using a wounded skin model.**

## Prototypes

Superabsorbent wound dressings (Prototype Dressings 1 based on Zetuvit plus<sup>®</sup>) and ointment impregnated wound contact layers (Prototype Dressings 2 based on Atrauman<sup>®</sup>) were coated with a Mixed Metal Complex (MMC). The metal complex is prepared by complexation of ionic Silver (Ag) and Zinc (Zn) with Ethylenediaminetetraacetic acid (EDTA). Dressing concentrations of 3% MMC for Prototype Dressing 1 and 1.8% MMC for Prototype Dressing 2 were used to design new antimicrobial wound dressing prototypes with antibiofilm properties against deep wounds infected with biofilm.

## Methodology


Intact human skin tissue (HST) explants (5 mm) were produced with 2 mm partial thickness wounds by biopsy punch. Explant wound beds were inoculated with *Pseudomonas aeruginosa* or methicillin-resistant *Staphylococcus aureus* (Xen30) and incubated to establish 48 or 72 hour biofilms.



Wounds beds were washed to remove planktonic bacteria and treated with liquid 3% MMC active solution (100  $\mu$ L/treatment), commercially available Zetuvit plus (Dressing 1)\* and non-medicated wound contact layer Atrauman (Dressing 2)#, or respective new prototype dressings (Prototype Dressing 1 or Prototype Dressing 2).

After 24 hours of treatment, explants were sonicated, vortexed, and released bacteria plated for enumeration. Biofilm viability was assessed by Live/Dead staining of one explant replicate and imaged with an Olympus BX63 fluorescent microscope. Two explants from each treatment were reserved for histological evaluation with hematoxylin and eosin (H&E) and Gram-Twort staining.


## Results




Organism	Log Reduction from Growth Control (Mean $\pm$ SD)				
	Liquid Active	Dressing 1*	Prototype Dressing 1	Dressing 2#	Prototype Dressing 2
<i>P. aeruginosa</i> (ATCC <sup>®</sup> 700888 <sup>™</sup> )	>3.35 $\pm$ 0.30	0.36 $\pm$ 0.05	3.56 $\pm$ 0.42	0.39 $\pm$ 0.08	1.46 $\pm$ 0.07
MRSA (Xen30)	>2.98 $\pm$ 0.96	0.67 $\pm$ 0.02	0.98 $\pm$ 0.25	0.30 $\pm$ 0.02	0.62 $\pm$ 0.08

**Table 1:** Average log reduction (CFU/explant) from growth control (n=3, Mean $\pm$ SD). *P. aeruginosa* 72h growth control = 8.28  $\pm$  0.01 CFU/explant. MRSA 96h growth control = 7.41  $\pm$  0.07 CFU/explant. \*Zetuvit<sup>®</sup> Plus, Hartmann, 2023. # Atrauman<sup>®</sup>, Hartmann, 2023.

Newly formulated Prototype Dressing 1 demonstrated the greatest overall percent reduction for both organisms, 99.97% reduction from *P. aeruginosa* 48 hour controls and 89.47% reduction from MRSA 96 hour controls.



**Figure 1:** Live/Dead fluorescent images of wounded human skin tissue. A) Uninfected HST wound bed after 72 hours of incubation. B) Wound bed 72 hours post-infection with *P. aeruginosa* ATCC<sup>®</sup> 700888<sup>™</sup> (end of treatment period). C) Wound bed 96 hours post-infection with *S. aureus* Xen 30 (end of treatment period). (Green: SYTO 9 stain, alive cells; red: propidium iodide, dead cells).



**Figure 2:** Haematoxylin & Eosin and Gram-Twort staining images of wounded HST infected with MRSA Xen 30, with and without treatment.

## Discussion and Conclusions

- Greater efficacy for Prototype Dressing 1 compared to legacy Dressing 1\* was supported by log CFU/explant recovered post treatment as well as visibly less biofilm remaining in wound beds of tissue explants as observed by histology. Prototype Dressing 2 was likewise more effective than legacy impregnated Dressing 2#.
- Newly designed prototype dressings demonstrated greater antibiofilm efficacy than their existing commercial products.** Both prototype dressings were outperformed only by the liquid active solution alone.
- H&E and Gram-Twort stained histology showed that the bacteria, although killed, remained in the wound bed when treated with liquid active alone as compared to any of the physical dressings assessed.
- Physical wound dressings may function to physically remove bacteria from the wound bed when discarded post treatment.